



Axon Spectrashield Tinted Soft Contact Lenses

Lens Information

Lens Requested: Axon Spectrashield Tinted Soft Contact Lens

Lens Material: Hioxifilcon 49%

Lens Source: Tint is applied to Special Eyes brand Hioxifilcon 49% lens.

Lens Type: Spherical and Toric Rx only

Lens Duration: 365 day lens

Lens Parameters

Diameter: 13.5 – 16.0 mm in 0.1 steps

Base Curve: 7.8 - 9.2 in 0.1 steps

Sphere Power: -25.0D to +25.0D in 0.1D steps

Toric Power: 15.0D

Cylinder Power Range: 0° to 360° in 1° steps

Tint: FDA Cleared. 510(k) number K991995



CONTACT LENSES ARE MADE TO ORDER BASED ON
PATIENT'S PRESCRIPTION

